



Vannas-Tübingen Spring Scissors - Angled Up/Sharp/8.5cm/5mm

Cutting Edge



Allgemein

Produktnummer	15005-08
Kategorie	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504908



Finde 15005-08 auf
finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	3 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	8.5 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		