



Vannas-Tübingen Spring Scissors - Angled Up/Sharp/8.5cm/5mm
Cutting Edge



General Info

Product Number	15005-08
Category	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504908



Find 15005-08 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	8.5 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		