



Castroviejo Spring Scissors - Angled to Side/Sharp/10cm/10mm Cutting Edge



General Info

|                |                                                            |
|----------------|------------------------------------------------------------|
| Product Number | 15006-09                                                   |
| Category       | Scissors & Micro Scissors > Micro Scissors<br>MedTech Line |
| EAN            | 4251909504915                                              |



Find 15006-09 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |     |
|------------------------|---------------------------------|-----------------|-----|
| Shipping Weight        | 30 g                            | Net Weight      | 5 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA  |
| Minimum Order Quantity | 1                               |                 |     |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |     |
| Cuttingedge            | 10 mm                           |                 |     |
| Length                 | 10 cm                           |                 |     |
| Serrated               | No                              |                 |     |
| Tipdiameter            | 0.1 mm                          |                 |     |
| Tips                   | Sharp/Sharp                     |                 |     |
| Tipshape               | Angled to side                  |                 |     |