



Vannas Spring Scissors - Microserrated/Straight/Sharp/9cm/8mm

Cutting Edge



Microserrated to minimize slippage.



General Info

Product Number	15007-08
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909504922



Find 15007-08 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	8 mm		
Length	9 cm		
Serrated	Yes		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		