



Vannas Spring Scissors - Straight/Sharp/9cm/8mm Cutting Edge



Ultra long cutting edge without sacrificing tip width.

General Info

Product Number	15009-08
Category	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909504946



Find 15009-08 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	8 mm		
Length	9 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		