



Spring Scissors - Angled up/10cm/4mm Cutting Edge



These spring scissors aid in cutting tissues from various angles.

General Info

Product Number	15010-10
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909504960



Find 15010-10 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	4 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		