



Noyes Spring Scissors - Angled Up/Sharp/12cm/14mm Cutting Edge



Allgemein

Produktnummer	15013-12
Kategorie	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909504991



Finde 15013-12 auf
finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	11 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	14 mm		
Length	12 cm		
Serrated	No		
Tipdiameter	0.2 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		