



Noyes Spring Scissors - Angled Up/Sharp/12cm/14mm Cutting Edge



General Info

Product Number	15013-12
Category	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909504991



Find 15013-12 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	14 mm		
Length	12 cm		
Serrated	No		
Tipdiameter	0.2 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		