



Wescott Spring Scissors - Slightly Curved Up/Sharp/11cm/15mm Cutting Edge



Allgemein

|               |                                            |
|---------------|--------------------------------------------|
| Produktnummer | 15015-11                                   |
| Kategorie     | Scissors & Micro Scissors > Micro Scissors |
| EAN           | 4251909505004                              |



Finde 15015-11 auf [finescience.com](https://finescience.com)

Technische Informationen

|                     |                                 |              |      |
|---------------------|---------------------------------|--------------|------|
| Versandgewicht      | 35 g                            | Nettogewicht | 11 g |
| Paketgröße          | 1 Pc.                           | Maßeinheit   | EA   |
| Mindestbestellmenge | 1                               |              |      |
| Alloymaterial       | Surgical Stainless Steel (Inox) |              |      |
| Cuttingedge         | 15 mm                           |              |      |
| Length              | 11 cm                           |              |      |
| Serrated            | No                              |              |      |
| Tipdiameter         | 0.2 mm                          |              |      |
| Tips                | Sharp/Sharp                     |              |      |
| Tipshape            | Curved                          |              |      |