

F · S · T

## Castroviejo Spring Scissors - Curved/Sharp/11cm/10mm Cutting Edge



## Allgemein

Produktnummer	15017-10
Kategorie	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909505011



Finde 15017-10 auf  
[finescience.com](https://www.finescience.com)

## Technische Informationen

Versandgewicht	30 g	Nettogewicht	10 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Cuttingedge	10 mm		
Length	11 cm		
Serrated	No		
Tipdiameter	0.2 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		