



Castroviejo Spring Scissors - Curved/Sharp/11cm/10mm Cutting Edge



General Info

Product Number	15017-10
Category	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909505011



Find 15017-10 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	10 mm		
Length	11 cm		
Serrated	No		
Tipdiameter	0.2 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		