




Vannas Spring Scissors - Straight/Sharp/10cm/4mm Cutting Edge



General Info

Product Number	15018-10
Category	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909505028



Find 15018-10 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	4 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		