



Vannas Spring Scissors - Curved/Sharp/10cm/4mm Cutting Edge



Allgemein

Produktnummer	15019-10
Kategorie	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909505035



Finde 15019-10 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	10 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	4 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		