



Spring Scissors - Straight/Sharp/14cm/6mm Cutting Edge



Allgemein

Produktnummer	15020-15
Kategorie	Scissors & Micro Scissors > Micro Scissors
EAN	4251826395511



Finde 15020-15 auf finescience.com

Technische Informationen

Versandgewicht	50 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	6 mm		
Length	14 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		