



Spring Scissors - Curved/Sharp/10cm/7mm Cutting Edge



Allgemein

Produktnummer	15023-10
Kategorie	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909505042



Finde 15023-10 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	7 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		