



Spring Scissors - Straight/Blunt/10cm/7mm Cutting Edge



Allgemein

Produktnummer	15025-10
Kategorie	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909505066



Finde 15025-10 auf
finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	7 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.3 mm		
Tips	Blunt/Blunt		
Tipshape	Straight		