



Spring Scissors - Straight/Blunt/10cm/7mm Cutting Edge



General Info

Product Number	15025-10
Category	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909505066



Find 15025-10 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.3 mm		
Tips	Blunt/Blunt		
Tipshape	Straight		