



Micro Scissors - Straight/Sharp/14.5cm/2.5mm Cutting Edge



General Info

Product Number	15030-14
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505097



Find 15030-14 on finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	2.5 mm		
Length	14.5 cm		
Serrated	Yes		
Tipdiameter	0.3 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		