



Biemer Vessel Spring Scissors - Straight/Blunt/8.5cm/4mm Cutting Edge



These scissors are designed to securely cut vessels without slipping.

General Info

Product Number	15032-08
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505110



Find 15032-08 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	4 mm		
Length	8.5 cm		
Serrated	No		
Tipdiameter	2 mm		
Tips	Blunt/Blunt		
Tipshape	Straight		