



Biemer Vessel Spring Scissors - Straight/Blunt/13cm/2mm Cutting Edge



These scissors are designed to securely cut vessels without slipping.

This product is discontinued and has limited availability.

General Info

Product Number	15032-13
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505127



Find 15032-13 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	2 mm		
Length	13 cm		
Serrated	No		
Tipdiameter	2 mm		
Tips	Blunt/Blunt		
Tipshape	Straight		