



Spring Scissors - Angled to Side/Sharp-Ball Tip/10cm/8mm Cutting Edge



The ball tip prevents piercing of underlying tissue.

General Info

Product Number	15033-09
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505134



Find 15033-09 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	8 mm		
Length	10 cm		
Feature	1.3 mm Ball Diameter		
Serrated	No		
Tipdiameter	1.3 mm		
Tips	Sharp/Ball tip		
Tipshape	Angled to side		