

F · S · T

Micro Scissors - Straight/Sharp/14cm/2.2mm Cutting Edge

For the finest microdissection, these spring scissors have the smallest blades. The black handles are made of ribbed plastic to ensure a solid grip.

General Info

Product Number	15036-14
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505141



Find 15036-14 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	55 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Cuttingedge	2.2 mm		
Length	14 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		