



Dowell Scissors - Angled Up/11.5cm/10mm Cutting Edge



Allgemein

Produktnummer	15042-11
Kategorie	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505172



Finde 15042-11 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	13 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	10 mm		
Length	11.5 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		