



Dowell Scissors - Angled Up/11.5cm/10mm Cutting Edge



General Info

Product Number	15042-11
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505172



Find 15042-11 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	10 mm		
Length	11.5 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Angled up		