



Micro 2000 Spring Scissors - ToughCut®/Curved/Sharp/12.5cm/6mm
Cutting Edge



Allgemein

Produktnummer	15123-12
Kategorie	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909505226



Finde 15123-12 auf
finescience.com

Technische Informationen

Versandgewicht	50 g	Nettogewicht	13 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Cuttingedge	6 mm		
Length	12.5 cm		
Feature	ToughCut®		
Serrated	Yes		
Tipdiameter	0.3 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		