



Micro 2000 Spring Scissors - ToughCut®/Straight/Sharp/12.5cm/6mm
Cutting Edge



General Info

Product Number	15124-12
Category	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909505233



Find 15124-12 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	6 mm		
Length	12.5 cm		
Feature	ToughCut®		
Serrated	Yes		
Tipdiameter	0.3 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		