



Large Moria Vannas Spring Scissors - Straight/Sharp/10cm/7mm Cutting Edge



Allgemein

Produktnummer	15370-52
Kategorie	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505264



Finde 15370-52 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	9 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.15 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		