



Large Moria Vannas Spring Scissors - Straight/Sharp/10cm/7mm Cutting Edge



General Info

Product Number	15370-52
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505264



Find 15370-52 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.15 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		