



Moria Spring Scissors - Curved/Sharp/8cm/5mm Cutting Edge



These have the finest tip diameter of all the Moria scissors we sell. Hand fashioned using 40x magnification, these tips are especially useful when incising vessel walls for catheter insertion.

General Info

Product Number	15396-01
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505295



Find 15396-01 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	8 cm		
Serrated	No		
Tipdiameter	0.15 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		