



Vannas Round Handled Spring Scissors - Curved/Sharp/12.5cm/6mm
Cutting Edge



Round handles allow easy fingertip adjustments.

General Info

Product Number	15401-12
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505318



Find 15401-12 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	6 mm		
Length	12.5 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		