



Vannas Spring Scissors - Titanium/Straight/Sharp/9.5cm/5mm Cutting Edge



Allgemein

Produktnummer	15610-08
Kategorie	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505332



Finde 15610-08 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	3 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Cuttingedge	5 mm		
Length	9.5 cm		
Serrated	No		
Tipdiameter	0.3 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		