



Vannas Spring Scissors - Titanium/Curved/Sharp/9.5cm/5mm Cutting Edge



General Info

Product Number	15615-08
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505349



Find 15615-08 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	3 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Cuttingedge	5 mm		
Length	9.5 cm		
Serrated	No		
Tipdiameter	0.3 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		