



Ceramic Coated Titanium - Curved/11cm/5mm Cutting Edge



General Info

Product Number	15651-11
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909505370



Find 15651-11 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	50 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Coating	Ceramic		
Antimagnetic	Yes		
Cuttingedge	5 mm		
Length	11 cm		
Serrated	No		
Tipdiameter	0.2 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		