



Ceramic Coated Spring Scissors - Straight/10cm/7mm Cutting Edge



Allgemein

Produktnummer	15750-11
Kategorie	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909505387



Finde 15750-11 auf
finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	7 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Ceramic coated Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		