



Ceramic Coated Spring Scissors - Straight/10cm/7mm Cutting Edge



General Info

Product Number	15750-11
Category	Scissors & Micro Scissors > Micro Scissors MedTech Line
EAN	4251909505387



Find 15750-11 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Ceramic coated Surgical Stainless Steel (Inox)		
Cuttingedge	7 mm		
Length	10 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		