



Vannas Spring Scissors - Highly Polished/Straight/Sharp/8.5cm/5mm
Cutting Edge



The highly polished surface makes the scissors more resistant to corrosion and easier to clean. It also prevents tissue from sticking on it while operating, it's also called non-sticky surface.

General Info

Product Number	15903-08
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909545215



Find 15903-08 on
finescience.com

Technical Info

Shipping Weight	24 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	8.5 cm		
Feature	High polished		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		