



Mini-Friedmann Rongeur - Curved/2.5mm Cup Diameter/12.5cm



Strong and large jaws are ideal for use on the harder, outer regions of the spinal cord as well as for gaining access to delicate structures encased in bone.

General Info

|                |                                         |
|----------------|-----------------------------------------|
| Product Number | 16000-14                                |
| Category       | Rongeurs, Punches & Sampling > Rongeurs |
| EAN            | 4251909505400                           |



Find 16000-14 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 70 g                            | Net Weight      | 70 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Cupsize                | 2.5 mm                          |                 |      |
| Length                 | 12.5 cm                         |                 |      |
| Jointaction            | Single                          |                 |      |
| Tipshape               | Curved                          |                 |      |