



Luer Rongeur - Curved/4mm Cup Diameter/15cm



General Info

Product Number	16001-15
Category	Rongeurs, Punches & Sampling > Rongeurs
EAN	4251909505417



Find 16001-15 on
finescience.com

Technical Info

Shipping Weight	110 g	Net Weight	107 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cupsize	4 mm		
Length	15 cm		
Jointaction	Single		
Tipshape	Curved		