



Beyer Rongeur - Curved/2.5mm Cup Diameter/18cm



The double action joint relieves strain on your hand.

General Info

Product Number	16002-18
Category	Rongeurs, Punches & Sampling > Rongeurs
EAN	4251909505424



Find 16002-18 on
finescience.com

Technical Info

Shipping Weight	135 g	Net Weight	127 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cupsize	2.5 mm		
Length	18 cm		
Jointaction	Double		
Tipshape	Curved		