



Pearson Rongeur - Shallow Bend/1mm Cup Diameter/17cm



Provides excellent visibility at dissection site due to its extra long and narrow approach to the fine cutting jaws. Especially suitable for larger hands. Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.

General Info

Product Number	16015-17
Category	Rongeurs, Punches & Sampling > Rongeurs
EAN	4251909505431



Find 16015-17 on
finescience.com

Technical Info

Shipping Weight	110 g	Net Weight	103 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Cupsize	1 mm		
Length	17 cm		
Jointaction	Single		
Tipshape	Shallow bend		