



Pearson Rongeur - Sharp Bend/1.2mm Cup Diameter/17cm



Provides excellent visibility at dissection site due to its extra long and narrow approach to the fine cutting jaws. Especially suitable for larger hands. Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.

General Info

|                |                                         |
|----------------|-----------------------------------------|
| Product Number | 16016-17                                |
| Category       | Rongeurs, Punches & Sampling > Rongeurs |
| EAN            | 4251909505448                           |



Find 16016-17 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |       |
|------------------------|---------------------------------|-----------------|-------|
| Shipping Weight        | 105 g                           | Net Weight      | 104 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA    |
| Minimum Order Quantity | 1                               |                 |       |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |       |
| Cupsize                | 1.2 mm                          |                 |       |
| Length                 | 17 cm                           |                 |       |
| Jointaction            | Single                          |                 |       |
| Tipshape               | Sharp bend                      |                 |       |