



Kleinert-Kutz Rongeur - Curved/2mm Cup Diameter/15cm



The double action joint relieves strain on your hand.

General Info

Product Number	16022-15
Category	Rongeurs, Punches & Sampling > Rongeurs
EAN	4251909505479



Find 16022-15 on
finescience.com

Technical Info

Shipping Weight	95 g	Net Weight	91 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cupsize	2 mm		
Length	15 cm		
Jointaction	Double		
Tipshape	Curved		