



Kleinert-Kutz Rongeur - Straight/2mm Cup Diameter/15cm



The double action joint relieves strain on your hand.

General Info

Product Number	16023-15
Category	Rongeurs, Punches & Sampling > Rongeurs
EAN	4251909505486



Find 16023-15 on
finescience.com

Technical Info

Shipping Weight	100 g	Net Weight	98 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cupsize	2 mm		
Length	15 cm		
Jointaction	Double		
Tipshape	Straight		