




Bone Pliers - Straight/15.5cm



With smooth, flat tips. These are ideal for peeling back skull flaps.

General Info

Product Number	16025-14
Category	Clamps & Clips
EAN	4251909505493



Find 16025-14 on
finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	62 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Jointaction	Single		
Tips	Flat		
Tipshape	Straight		