



Bone Scissors - Tungsten Carbide/Straight/12.5cm/15mm Cutting Edge



Strong heavy blades for cutting bone.



General Info

Product Number	16030-12
Category	Scissors & Micro Scissors > Surgical Scissors
EAN	4251909505509



Find 16030-12 on finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	45 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Cuttingedge	15 mm		
Length	12.5 cm		
Feature	Tungsten Carbide		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		