



Bone Scissors - Tungsten Carbide/Straight/16cm/15mm Cutting Edge



Strong heavy blades for cutting bone.



General Info

Product Number	16030-16
Category	Scissors & Micro Scissors > Surgical Scissors
EAN	4251909505516



Find 16030-16 on finescience.com

Technical Info

Shipping Weight	75 g	Net Weight	70 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	15 mm		
Length	16 cm		
Feature	Tungsten Carbide		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		