



Skull Flap Cutter - Angled Up/14cm/10mm Cutting Edge



The 15 angled up tips have a narrow and thin platform to fit between dura and skull bone. The double action joint relieves strain on your hand.

General Info

Product Number	16050-15
Category	Cutters > Medical
EAN	4251909505523



Find 16050-15 on finescience.com

Technical Info

Shipping Weight	90 g	Net Weight	88 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladematerial	Stainless Steel (Inox)		
Cuttingedge	10 mm		
Length	14 cm		
Jointaction	Double		
Tipshape	Angled up		