



Bone Forceps - Angled/Cross Serrated/11.5



Strong cross serrated tips for bone crushing and removal.

General Info

Product Number	16060-11
Category	Forceps & Tweezers > Serrated
EAN	4251909505530



Find 16060-11 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Tipdiameter	1 mm		
Tipdimensions	1 x 1 mm		
Tips	Cross serrated		
Tipshape	Angled		