



Friedmann-Pearson Rongeur - Straight/0.7mm Cup Diameter/14cm



Best suited for delicate work, these are the finest tipped rongeurs in our product line. The fine-cupped jaws have been designed with a longer and slimmer approach for better visibility. Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.

General Info

Product Number	16120-14
Category	Rongeurs, Punches & Sampling > Rongeurs Most Popular > Rongeurs, Punches & Sampling
EAN	4251909505639



Find 16120-14 on
finescience.com

Technical Info

Shipping Weight	60 g	Net Weight	52 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cupsize	0.7 mm		
Length	14 cm		
Jointaction	Single		
Tipshape	Straight		