



Lempert Cutter - Straight/11cm/9mm Cutting Edge



One edge is finely serrated.



General Info

Product Number	16140-11
Category	Cutters > Medical
EAN	4251909505660



Find 16140-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	9 mm		
Length	11 cm		
Jointaction	Single		
Serrated	Yes		
Tipshape	Straight		